

# Perfect LL Gel System W

■ Perfect LL Gel W (カセットサイズ:160W×100Hmm / ゲルサイズ:140W×80H×1.0Tmm) 保存:4°C(凍結不可)  
タンパク質分析用 (Tris-HCL 中性ゲル) SDS-PAGE / Native PAGE

製品名		カタログ番号	製品内容		標準価格
ゲルサイズ	ゲルバッファー		ゲル濃度 (%)	well	10枚入り
Perfect LL Gel W タンパク質分析用	Tris-HCL 中性	DLL-706P	5%	20well	¥33,300
		DLL-707P	5%	28well	¥33,300
		DLL-702P	5%	2-D	¥33,300
		DLL-716P	7.5%	20well	¥33,300
		DLL-717P	7.5%	28well	¥33,300
		DLL-712P	7.5%	2-D	¥33,300
		DLL-726P	10%	20well	¥33,300
		DLL-727P	10%	28well	¥33,300
		DLL-722P	10%	2-D	¥33,300
		DLL-736P	12.5%	20well	¥33,300
		DLL-737P	12.5%	28well	¥33,300
		DLL-732P	12.5%	2-D	¥33,300
		DLL-746P	15%	20well	¥33,300
		DLL-747P	15%	28well	¥33,300
		DLL-742P	15%	2-D	¥33,300
		DLL-756P	5-10%	20well	¥33,300
		DLL-757P	5-10%	28well	¥33,300
		DLL-752P	5-10%	2-D	¥33,300
		DLL-7E6P	5-12.5%	20well	¥33,300
		DLL-7E7P	5-12.5%	28well	¥33,300
		DLL-7E2P	5-12.5%	2-D	¥33,300
		DLL-766P	5-15%	20well	¥33,300
		DLL-767P	5-15%	28well	¥33,300
		DLL-762P	5-15%	2-D	¥33,300
		DLL-776P	5-20%	20well	¥33,300
		DLL-777P	5-20%	28well	¥33,300
		DLL-772P	5-20%	2-D	¥33,300
		DLL-786P	7.5-15%	20well	¥33,300
		DLL-787P	7.5-15%	28well	¥33,300
		DLL-782P	7.5-15%	2-D	¥33,300
DLL-796P	10-20%	20well	¥33,300		
DLL-797P	10-20%	28well	¥33,300		
DLL-792P	10-20%	2-D	¥33,300		
DLL-7X6P	特注濃度	20well	¥33,300		
DLL-7X7P	特注濃度	28well	¥33,300		
DLL-7X2P	特注濃度	2-D	¥33,300		

## ■ Perfect LL Gel W 用 装置・オプション

カタログ番号	製品名	入数	標準価格
PCE-00WC	ERICA-W 電気泳動槽(ダミープレート付)	1台	¥168,000
PCE-PLEW	ERICA-W 用ダミープレート	1枚	¥3,000
DRC-W6LG	Sample loading Guide 20well 用	1枚	¥4,200
DRC-W7LG	Sample loading Guide 28well 用	1枚	¥4,200